Partition panels

SURGICAL - PASS-THROUGH CABINET FOR OPERATING



Pass-through cabinets are used to supply materials to the operating theatre from the sterile hallway. The cabinets are produced in two versions. One version is equipped with three stainless steel shelves and with a travel-through space at the bottom. The travel-through space allows personnel to transport stainless steel containers with sterile materials (or sterile materials without containers) from the central sterilisation room on medical carts of standardized dimensions. One person opens the pass-through cabinet from the sterile hallway and pushes the cart inside. Another person from the operating theatre then takes the cart out. Any remaining material can be stored on shelves in the pass-through cabinet. The second version is equipped with six shelves. It is recommended that each operating theatre have one or two pass-through cabinets in order to create sufficient storage space for materials necessary for the surgical procedures.

TECHNICAL DATA

Pass-through cabinet for operating theatres
Pass-through cabinet

Dimensions			
width	depth	height	
1020 mm	830 mm	2100 mm	

Version		
Cabinet with shelves for operating theatre		
Travel-through cabinet for operating theatre		

Blocker

Cabinet with indication

Cabinet with indication and blocker

Without indication

Glazing

Single-sided antireflective glazing

Double-sided antireflective glazing

Single-sided safety glazing

Double-sided safety glazing

Without glazing

MORE INFORMATION, PHOTOS

